



**SNS COLLEGE OF ALLIED HEALTH SCIENCES**  
SNS Kalvi Nagar, Coimbatore - 35  
Affiliated to Dr MGR Medical University, Chennai



**DEPARTMENT OF CARDIO PULMONARY PERFUSION CARE**  
**TECHNOLOGY**

**COURSE NAME : PRINCIPLES OF PERFUSION**

**2<sup>ND</sup> YEAR**

**TOPIC : ASSEMBLING OF THE CIRCUIT**



# Safety Checks and Setup of the Circuit



## Checking the heart–lung machine and accessories

### 1: Connection checks:

- Checking the cables, plugs, and sockets.
- All cables should be placed away from movement to avoid any accidents, or cable damage.
- Checking the power supply for the heater-cooler and pump light.
- Connect the gas hoses to the gas source, and Checking the gas supplies of air and oxygen
- Checking the mixers flow meters for leaks.








# Safety Checks and Setup of the Circuit (cont)



## 2: Pump head checks:

- Power on.
- The pump heads are free from foreign bodies.
- The Rollers and guides move smoothly.
- The pump heads rotate in the correct direction.
- Calibrate flow scale rpm to lpm on control panels for roller pumps.





## 2: Pump head checks (cont)



- Availability of winding handles in case of power failure is a must.
- The size of tubing inserts are correct for the tubes size used.
- The tubing should be aligned within the pump housing.
- Properly assembled roller pump tubing must resembles a "**U**" **shape**, advantages of that are: avoiding more tight non needed occlusion, lowering the load on head pump, decreasing spallation and hemolysis






## 3: Electrical safety device checks:



- Battery backup (UPS) is charged.
- Turn on the pump heads, and unplug the power cable of the CPB machine to ascertain that the power failure alarm, and backup battery are working properly.
- Pressure transducers are calibrated.
- Level detectors are working properly .
- Bubble detectors are working properly





# Setup of disposable heart lung equipment



## Set up time

- The responsible Perfusionist takes **10-15 minutes** to set up a CPB circuit.

## Set up shelf life

- The set up can be kept on standby for up to **3 weeks without priming** (Dry).
- **When assembling CPB circuits using a sterile technique, all ports should be sealed (closed) tightly with provided caps, and water lines are not connected with the oxygenator.**



# Setup of disposable heart lung equipment (cont)



## 1. The Hard shell Venous Reservoir setup:

- A. It should be check expire dates.
- B. It should be check sterility for pack (if the pack opened or damaged don't use the materials inside), then removed it from the package and examined for any defects.
- C. Place the venous reservoir on its holder and orient it to allow full view of the reservoir
- D. Remove the yellow gas vent port cap located on the top of the reservoir.
- E. Connect the pump header line to the venous reservoir outlet and oxygenator (membrane) inlet.
- F. Connect the venous line to the venous inlet port.
- G. Connect any suction port required. (**Warning: Unused ports should be sealed (closed) tightly with provided caps.**)
- H. Connect the recirculation line to the recirculation port on the venous inlet.





## 2: The oxygenator (membranes) setup:



- A. It should be check expire dates.
- B. It should be check sterility for pack (if the pack opened or damaged don't use the materials inside), then removed it from the package and examined for any defects.
- C. Place the membrane carefully on its holder. Be sure that the venous reservoir is higher than the membrane.
- D. Connect the water inlet and outlet lines to the heat exchanger water ports. Start water flow and check for leaks from the water compartment.
- E. Connect the gas line to the gas inlet port.
- F. Remove the cap from the gas outlet port.
- G. In integrated arterial filters without self-venting the premembrane air purge line must remain open to ensure proper air removal.





### 3: Arterial line filters setup if demanded:



- A. It should be check expire dates.
  - B. It should be check sterility for pack (if the pack opened or damaged don't use the materials inside),then removed it from the package and examined for any defects.
  - C. Place the arterial filter carefully on its holder avoiding any kinking of lines.
  - D. Connect the pressure line to the vent port of the arterial filter (on 3-ways stopcock) and connect the end of the pressure line to the pressure transducer inlet.
  - E. Enter the pressure setting (warning and limit) in the control monitor.
- 
- ❖ If the pressure exceeds a **set warning** (usually 300 mmHg) **the pump will alarming.**
  - ❖ If the line pressure exceeds **a set limit** (usually 350 mmHg), **the pump will stop.**



## 4: The cardioplegia system setup if demanded:



- A. It should be check expire dates.
- B. It should be check sterility for pack (if the pack opened or damaged don't use the materials inside), then removed it from the package and examined for any defects.
- C. Place the cardioplegia carefully on its holder.
- D. Connect the water lines to the cardioplegia administration set, and start water flow to ensure that it is free from leaks.
- E. Assemble circuit according to the instructions from the manufacturer's manual.
- F. The cardioplegia pressure transducer and purge lines are connected to the cardioplegia delivery device.



## 5: The centrifugal pump setup if demanded:



- A. It should be check expire dates.
- B. It should be check sterility for pack (if the pack opened or damaged don't use the materials inside), then removed it from the package and examined for any defects.
- C. Checked the head motor is free from dirt or foreign bodies to maintain the function of the device.
- D. Place the centrifugal cones carefully on the centrifugal head motor.
- E. Connect the PVC tube to the outlet of the venous reservoir and connect the other end of the tube to the inlet port of the centrifugal cones.
- F. Connect another tube to the outlet port of the centrifugal cones and connect the other end of the tube to the oxygenator inlet port.
- G. Ensure that the relevant handles-crank is available in case of power failure.





THANK YOU