

## Chapter D – Faradic Current

1. Fill in the Blank: Faradic current has a pulse duration of \_\_\_\_\_ ms.
2. Riddle: 'Interrupted, asymmetrical, I contract weak muscles.'
3. Odd-One-Out: Faradic, IDC, TENS, Ultrasound.
4. Jumbled Word: FDRCAIAA
5. Crossword Clue: Current type used for innervated weak muscle (7 letters).
6. True/False: Faradic current can stimulate denervated muscle.
7. Match the Following: Pulse duration → 0.1–1 ms, Frequency → 50 Hz, Type → Asymmetrical, Effect → Muscle contraction
8. Riddle: 'I'm used in foot drop re-education.'
9. Fill in the Blank: Faradic current is used for \_\_\_\_\_ stimulation (sensory/motor).
10. Odd-One-Out: Foot drop, Denervated muscle, Weak muscle, Facial palsy.